

ME 08 - BELL CRANK LEVER APPARATUS

EM&E



Apparatus consists of a wooden right angled lever with fulcrum at the angle, the horizontal arm is of 75cm and vertical arm is of 20cm. A spring balance with adjustable wing nut is attached to the shorter arm. The lever arm with scale is engraved at every 5cm. Apparatus is supplied with one sliding weight of 1 kg.

EDUCATIONAL SYSTEM YOU CAN BUY FROM:

ENGINEERING MODELS & EQUIPMENT

**Regd. Office:**

Mahavir Jain Market  
7 - Civil Lines, P.B. No. - 13  
Roorkee - 247667

Email: [emod@vsnl.com](mailto:emod@vsnl.com), [enggmodels@bsnl.in](mailto:enggmodels@bsnl.in)

**Factory & Works:**

Adarsh Nagar  
New Hardwar Road  
Roorkee - 247667 (U.K.)  
Phones: 01332-273721  
Website: [enggmod.com](http://enggmod.com)

AN ISO 9001:2000 COMPANY