

ME 09 - BENDING MOMENT APPARATUS



Apparatus consists of a wooden beam 180cm long hinged at a distance of 60cm from one end. A spring balance with adjusting screw compensates the initial extension of the spring and keeps the beam horizontal under any load. Apparatus is supplied complete with three 'U' shaped sliding weights and two knife edge end supports on stands.

EDUCATIONAL SYSTEM YOU CAN BUY FROM:

ENGINEERING MODELS & EQUIPMENT

Regd. Office:

Mahavir Jain Market
7 - Civil Lines, P.B. No. - 13
Roorkee - 247667

Email: emod@vsnl.com, enggmodels@bsnl.in

Factory & Works:

Adarsh Nagar
New Hardwar Road
Roorkee - 247667 (U.K.)
Phones: 01332-273721
Website: enggmod.com

AN ISO 9001:2000 COMPANY