

ME 13 - COMBINED COIL AND BELT FRICTION

EM&E



Apparatus consists of a metallic pulley of 30 cm diameter with two grooves, one for rope and the other for belt. At back side of the pulley, a circular wooden disc is provided with graduation at  $10^{\circ}$  intervals. A rotating pulley is fixed at the circumference of the scale. Apparatus is supplied complete with one belt and rope. The apparatus is mounted on a heavy vertical stand.

EDUCATIONAL SYSTEM YOU CAN BUY FROM:

ENGINEERING MODELS & EQUIPMENT

**Regd. Office:**

Mahavir Jain Market  
7 - Civil Lines, P.B. No. - 13  
Roorkee - 247667

Email: [emod@vsnl.com](mailto:emod@vsnl.com), [enggmodels@bsnl.in](mailto:enggmodels@bsnl.in)

**Factory & Works:**

Adarsh Nagar  
New Hardwar Road  
Roorkee - 247667 (U.K.)  
Phones: 01332-273721  
Website: [enggmod.com](http://enggmod.com)

AN ISO 9001:2000 COMPANY