

ME 19 - DEFLECTION OF BEAM APPARATUS

EM&E



Apparatus consists of a M.S. beam with graduation scale. Supplied complete with two knife edge end supports, a hanger with knife edge and a pointer. Apparatus is supplied with a dial gauge of 25 mm travel with magnetic base and set of weights with hanger.

EDUCATIONAL SYSTEM YOU CAN BUY FROM:

ENGINEERING MODELS & EQUIPMENT

Regd. Office:

Mahavir Jain Market
7 - Civil Lines, P.B. No. - 13
Roorkee - 247667

Email: emod@vsnl.com, enggmodels@bsnl.in

Factory & Works:

Adarsh Nagar
New Hardwar Road
Roorkee - 247667 (U.K.)
Phones: 01332-273721
Website: enggmod.com

AN ISO 9001:2000 COMPANY