

## ME 50 - SIMPLE JIB CRANE APPARATUS



Apparatus consists of a tubular compression balance, pivoted about an axis fitted to the base. The tie chain has an adjustable angle bracket and fitted with a 10 kg extension balance. Apparatus is complete on a wooden base with 1.2m jib.

### OPTIONAL ACCESSORIES REQUIRED

Set of weight

*EDUCATIONAL SYSTEM YOU CAN BUY FROM:*

## ENGINEERING MODELS & EQUIPMENT

**Regd. Office:**

Mahavir Jain Market  
7 - Civil Lines, P.B. No. - 13  
Roorkee - 247667

Email: [emod@vsnl.com](mailto:emod@vsnl.com), [enggmodels@bsnl.in](mailto:enggmodels@bsnl.in)

**Factory & Works:**

Adarsh Nagar  
New Hardwar Road  
Roorkee - 247667 (U.K.)  
Phones: 01332-273721  
Website: [enggmod.com](http://enggmod.com)

*AN ISO 9001:2000 COMPANY*