

ME 31 - JOINTED ROOF TRUSS APPARATUS

EME



A metallic apparatus, about 1.2m in length whose one leg is fitted on an unmovable bracket and the other on rollers. Spring balances show tension on the rod and the compression on rafters. Apparatus is supplied complete with a dynamometer, chain and pan.

OPTIONAL ACCESSORIES REQUIRED

Set of weight

EDUCATIONAL SYSTEM YOU CAN BUY FROM:

ENGINEERING MODELS & EQUIPMENT

Regd. Office:

Mahavir Jain Market
7 - Civil Lines, P.B. No. - 13
Roorkee - 247667

Email: emod@vsnl.com, enggmodels@bsnl.in

Factory & Works:

Adarsh Nagar
New Hardwar Road
Roorkee - 247667 (U.K.)
Phones: 01332-273721
Website: enggmod.com

AN ISO 9001:2000 COMPANY