

ME 45 - SCREW JACK APPARATUS

ENG



Apparatus consists of a metallic machined screw with a pitch of 5mm carrying a double flanged turned table of about 20cm dia. Screw jack is fitted on a heavy C.I. base and is supplied complete with 2 adjustable pulleys, cord and hooks.

EDUCATIONAL SYSTEM YOU CAN BUY FROM:

ENGINEERING MODELS & EQUIPMENT

Regd. Office:

Mahavir Jain Market
7 - Civil Lines, P.B. No. - 13
Roorkee - 247667

Email: emod@vsnl.com, enggmodels@bsnl.in

Factory & Works:

Adarsh Nagar
New Hardwar Road
Roorkee - 247667 (U.K.)
Phones: 01332-273721
Website: engmod.com

AN ISO 9001:2000 COMPANY