


ME 48 - SHEAR FORCE APPARATUS

EME



Apparatus consists of two knife edge supports on stands and a bar 180cm long which is hinged at one point of its length. The dynamometer with an adjusting screw for restoring line of beam under load conditions is supplied with the apparatus. Apparatus is supplied complete with three sliding 'U' shaped weights.

EDUCATIONAL SYSTEM YOU CAN BUY FROM:

ENGINEERING MODELS & EQUIPMENT

Regd. Office:

Mahavir Jain Market
7 - Civil Lines, P.B. No. - 13
Roorkee - 247667

Email: emod@vsnl.com, enggmodels@bsnl.in

Factory & Works:

Adarsh Nagar
New Hardwar Road
Roorkee - 247667 (U.K.)
Phones: 01332-273721
Website: enggmod.com

AN ISO 9001:2000 COMPANY