


ME 59 - UNIVERSAL FORCE TABLE

EM&E



Apparatus consists of a circular aluminium disc of 40cm diameter, graduated in 360° . Apparatus is supplied complete with leveling screws, clamping device to fix the table at any desired angle, five sliding clamp pulleys, central ring, string and five sets of iron nickled slotted weights (each set containing nine weights and one hanger of 50 gm each).

EDUCATIONAL SYSTEM YOU CAN BUY FROM:

ENGINEERING MODELS & EQUIPMENT

Regd. Office:

Mahavir Jain Market
7 - Civil Lines, P.B. No. - 13
Roorkee - 247667

Email: emod@vsnl.com, enggmodels@bsnl.in

Factory & Works:

Adarsh Nagar
New Hardwar Road
Roorkee - 247667 (U.K.)
Phones: 01332-273721
Website: enggmod.com

AN ISO 9001:2000 COMPANY