

ME 61 - WINCH CRAB APPARATUS

EM&E



Apparatus consists of a heavy duty C.I. wall brackets. The grooved wheel is of 25 cm dia, and a set of gearing arrangements (machine cut gears) are provided for measurement of efficiency of winch crab.

OPTION

- (a) Single Purchase
- (b) Double Purchase

OPTIONAL ACCESSORIES REQUIRED

Set of weight

EDUCATIONAL SYSTEM YOU CAN BUY FROM:

ENGINEERING MODELS & EQUIPMENT

Regd. Office:

Mahavir Jain Market
7 - Civil Lines, P.B. No. - 13
Roorkee - 247667

Email: emod@vsnl.com, enggmodels@bsnl.in

Factory & Works:

Adarsh Nagar
New Hardwar Road
Roorkee - 247667 (U.K.)
Phones: 01332-273721
Website: enggmod.com

AN ISO 9001:2000 COMPANY