

**ME 15 - COMBINED INCLINED PLANE AND FRICTION SLIDE APPARATUS**

**ENGINEERING MODELS & EQUIPMENT**



Apparatus consists of a wooden plane 90 x 18 cm having adjustment for setting the angle precisely. Apparatus is supplied complete with a frictionless pulley, linear and circular scale, pan, a wheeled trolley and a set of eight slide draws.

**OPTIONAL ACCESSORIES REQUIRED**

Physical weight box

*EDUCATIONAL SYSTEM YOU CAN BUY FROM:*

**ENGINEERING MODELS & EQUIPMENT**

**Regd. Office:**

Mahavir Jain Market  
7 - Civil Lines, P.B. No. - 13  
Roorkee - 247667

Email: [emod@vsnl.com](mailto:emod@vsnl.com), [enggmodels@bsnl.in](mailto:enggmodels@bsnl.in)

**Factory & Works:**

Adarsh Nagar  
New Hardwar Road  
Roorkee - 247667 (U.K.)  
Phones: 01332-273721  
Website: [enggmod.com](http://enggmod.com)

*AN ISO 9001:2000 COMPANY*