

MODEL NO.
FMC-09

PELTON WHEEL TURBINE TEST RIG

ENGINEERING

EME



Technical Specification:

- ❖ Structure housing – Mild Steel.
- ❖ Sump tank (S.S.)
- ❖ Centrifugal pump – Three phase 5HP capacity.
- ❖ Control valves and drain valves.
- ❖ G.I pipe line

Experimental Capabilities:

- ❖ To determine the efficiency of the Pelton wheel turbine

Option:

Pelton Turbine Test Rig, 5 HP also available

EDUCATIONAL SYSTEM YOU CAN BUY FROM:

ENGINEERING MODELS & EQUIPMENT

Regd. Office:

Mahavir Jain Market
7 – Civil Lines, P.B. No. – 13
Roorkee – 247667 (U.K.)

Email: emod@vsnl.com, enggmodels@bsnl.in

Factory & Works:

Adarsh Nagar
New Hardwar Road
Roorkee – 247667 (U.K.)
Phones: 01332-273721
Website: enggmod.com

AN ISO 9001:2000 & 14001:2004 COMPANY