

MODEL NO.  
FMC-06

# HYDRAULIC RAM TEST RIG



## Technical Specification:

- ❖ Structure housing – Mild Steel.
- ❖ Sump tank (S.S.)
- ❖ Centrifugal monoblock pump – Single phase 1hp capacity.
- ❖ Control valves and drain valves of Gunmetal.
- ❖ G.I pipe line

## Experimental Capabilities:

- ❖ To study the Hydraulic ram.
- ❖ To determine the efficiency of the hydraulic ram.

*EDUCATIONAL SYSTEM YOU CAN BUY FROM:*

## ENGINEERING MODELS & EQUIPMENT

### Regd. Office:

Mahavir Jain Market  
7 – Civil Lines, P.B. No. – 13  
Roorkee – 247667 (U.K.)  
Fax: 01332-271621  
Email: [enggmod@nde.vsnl.net.in](mailto:enggmod@nde.vsnl.net.in), [info@enggmod.com](mailto:info@enggmod.com)

### Factory & Works:

Adarsh Nagar  
New Hardwar Road  
Roorkee – 247667 (U.K.)  
Phones: 01332-273721  
Website: [enggmod.com](http://enggmod.com)

*AN ISO 9001:2000 & 14001:2004 COMPANY*