

MODEL NO.
FM-07

APPARATUS CONDUCTING FOR ORIFICE & MOUTHPIECE EXPERIMENT



Technical Specification:

- ❖ Structure housing – Mild Steel.
- ❖ Sump tank (S.S.)
- ❖ Collecting tank (S.S.)
- ❖ A supply tank of S.S. sheet with over flow arrangement and gauge glass tube is provided. The tank has the provision for fixing the various orifices and mouth pieces (inter changeable).
- ❖ To measure the X-Y coordinates of jet coming out of orifice/mouthpiece, scale and sliding apparatus.
- ❖ Centrifugal monoblock pump – Single phase 0.5hp capacity.
- ❖ Control valves and drain valves of Gunmetal.
- ❖ G.I pipe line
- ❖ Strong M.S stand.

NOTE: APPARATUS SUBJECTED TO BE CHANGE AS PER DEMAND

EDUCATIONAL SYSTEM YOU CAN BUY FROM:

ENGINEERING MODELS & EQUIPMENT

Regd. Office:

Mahavir Jain Market
7 – Civil Lines, P.B. No. – 13
Roorkee – 247667 (U.K.)
Fax: 01332-271621
Email: emod@vsnl.com, enggmodels@bsnl.in

Factory & Works:

Adarsh Nagar
New Hardwar Road
Roorkee – 247667 (U.K.)
Phones: 01332-273721
Website: enggmod.com

AN ISO 9001:2000 & 14001:2004 COMPANY